

**LIVE**  
 MEASUREMENTS

**NEW**  
 DESIGN

# MOTICAM 1080N INT

Full HD Standalone camera for BA Series

The Moticam 1080N INT, designed for the BA Series, outputs fast live images in Full HD (1920x1080px) and displays your live images on an integrated LCD screen by connecting the camera through the HDMI port. Make live measurements, notes, image adjustments and more, thanks to its on-board software that can be controlled through the screen by connecting a USB mouse to the camera. By having an integrated micro SD card slot, saving images won't require a computer, making it easy to save and organize your images and videos. The USB port allows you to connect with a computer, even simultaneously, and use the Motic Images Plus software.

